



Nova5000 v1.00.XX EX

**Release Notes** (updated Sep 1, 2006)

**Firmware version:** **1.00.08 EX**  
**Release date:** August 8, 2006

**New Features:**

- Ability for Nova5000 data logger to communicate with MultiLab PC version
- Windows CE QFE (June 2006)
- ActiveSync files sync option added
- Nova Paint added
- Flash 7 added
- MultiLab CE 1.07 with upgraded toolbar icons and bug fixes
- Arial and Verdana fonts added
- Wordpad and TextMaker - default font size changed to 14
- CalcCE – Improved GUI
- Ability to mute all sounds added

**Bug Fixes:**

- Time/Date loss on reset fixed
- Suspend delay and Suspend failure fixed

**Additional Changes:**

Solitaire and FreeCell games were removed

**Firmware version:** **1.00.06 EX**  
**Release date:** June 15, 2006

**New Features:**

Microsoft monthly Windows CE 5.0 patch applied. Follow link for Microsoft release notes:

[http://download.microsoft.com/download/5/2/d/52dd9f6c-acf9-4089-9594f76524023eb2/Windows%20CE%205.0\\_Update\\_060531\\_2006M05.htm](http://download.microsoft.com/download/5/2/d/52dd9f6c-acf9-4089-9594f76524023eb2/Windows%20CE%205.0_Update_060531_2006M05.htm)

**Bug Fixes:**

Sound amplifier delay removed. Speaker shuts down immediately after the sound file stops playing (five second delay removed).

**Firmware version:** **1.00.05 EX**  
**Release date:** June 13, 2006

**New Features**

- VPN (PPTP) support
- Support for Direct Draw
- Support for Format Message API

**Bug Fixes:**

- Critical Error 608 resolved – Communication error when connecting PC USB cable to Nova5000 and attempting ActiveSync connection
- .CSV files now associated with PlanMaker
- Improved recognition of AC adapter connection via system tray icon